

Magnetic Resonance Imaging sequence programming

Zigzag EPI with Sinusoidal Gradient

Nan Lan



Method—Blipped EPI

Drawback:

- noise: round 110-120dB

Reason:

- frequent and sharp switch of gradient field




Fig. 1. Blipped Phase Encoding EPI with Sinusoidal Gradient
Source: <https://mriquestions.com>

Method—Sin Zigzag EPI

Pros:

- quite

Reason:

- smooth gradient trajectory due to sin gradient




Fig. 2. Zigzag-Aligned-Projections Phase Encoding EPI with Sinusoidal Gradient [2]

Image Result of Sin Zigzag EPI




Fig. 3. Image Result of Zigzag EPI with Sinusoidal Gradient

K-space result of Sin Zigzag EPI




Fig. 4. K-space Result of Zigzag EPI with Sinusoidal Gradient

Midpoint correction comparison





Fig. 5. Zigzag EPI with Sinusoidal Gradient without midpoint correction





Fig. 6. Zigzag EPI with Sinusoidal Gradient and midpoint correction

Com. between normal and zigzag sin EPI




Fig. 7. Nnormal EPI without Sinusoidal Gradient




Fig. 6. Zigzag EPI with Sinusoidal Gradient and midpoint correction

Reference

[2] Patrick Liebig, Robin M. Heidemann, Bernhard Hensel & David A. Porter *Accelerated silent echo-planar imaging*. ELSEVIER 2019. 55: 81–85.

